

Lissamine Green 1.5 mg Ophthalmic Strips

User Guidelines

Consumable

Intended use

Sterile, single-use dye strips for conjunctival and corneal vital staining during ocular surface examination.

Intended users

For use by trained healthcare professionals (ophthalmologists, optometrists and eye clinic staff). Follow your facility's protocols and applicable regulations.

Read this manual completely before first use and keep it with the product. This document provides general operating, care and safety guidance. Where the product is supplied with manufacturer labeling or instructions, those take precedence. For questions, a replacement copy, or a digital copy of this manual, email sales@scimedstore.com.

1. Safety

Safety warnings and cautions

- For use by trained healthcare professionals. Follow your facility's protocols and applicable regulations.
- Single use only. Use one strip per patient and discard after use. Do not reuse.
- Do not use if the individual wrapper is damaged or open, or if the expiry date has passed.
- Moisten only with sterile saline. Avoid touching the dye-impregnated tip with fingers.
- Do not use on patients with known sensitivity to the dye.
- Avoid touching the cornea with the strip.

2. Package contents

- 15 boxes of 100 individually wrapped sterile strips (1,500 strips total)

Check contents on receipt. Report missing or damaged items to sales@scimedstore.com within 48 hours of delivery.

3. Specifications

Dye	Lissamine green
Amount per strip	1.5 mg
Box size	100 strips
Lot size	15 boxes (1,500 strips)
Sterility	Sterile, individually wrapped
Use	Single use

4. Operation

Applying the strip

- 1 Check the expiry date and that the wrapper is intact.
- 2 Open the wrapper at the handle end without touching the dye tip.
- 3 Moisten the dye tip with sterile saline.
- 4 Ask the patient to look up and gently pull down the lower lid.
- 5 Touch the moistened tip lightly to the lower palpebral conjunctiva or fornix.
- 6 Ask the patient to blink gently to distribute the dye.
- 7 Examine with the slit lamp under white light and grade staining according to your clinical protocol.
- 8 Discard the strip immediately after use.

5. Cleaning and disinfection

- No cleaning required; the strips are single-use.
- Wipe any excess dye from the patient's lid margin with a clean tissue.

6. Maintenance and storage

- Store in the original boxes in a cool, dry place away from direct sunlight.
- Keep wrappers sealed until use.
- Check the expiry date on the box and wrapper before use; do not use expired strips.
- Use older stock first.

7. Disposal

- Dispose of used strips as clinical waste according to facility protocol.
- Dispose of expired or damaged strips with clinical waste.
- Recycle outer cartons where local rules permit.

8. Warranty and support

Products are covered by SciMed's standard warranty (12 months from purchase unless a different period is stated on your invoice). Warranty terms: scimedstore.com/pages/warranty.

Support, spare parts and additional copies of this manual: sales@scimedstore.com. Please include your order number and a photo of the product or issue.

SciMed Scientific Equipment LLC · San Jose, California, USA · scimedstore.com